


WATER TESTING
SOLUTIONS

Marine Microbiological Testing



Water Testing Solutions Ltd Unit 11 Wickham Business Centre, Harwood Road, Littlehampton, BN17 7AU, UK
Tel: +44 (0)1903 714199 www.watertestingolutions.com

Microbiological Testing

Potable water is our most important nutrient, used for both drinking and cooking. The water being used for personnel hygiene and all types of cleaning requires the same high quality. It is therefore important to have enough water of satisfactory quality to cover all types of usage. Safe potable water secures good health, but water can also contain hazardous microbes that cause infectious disease or food poisoning.

WTS-P101368 Aerobic Bacteria Test Pack (20)

- 20 x 2ml sterile syringes
- 20 x alcohol disinfecting wipes
- 20 x dilution pots
- 20 x compact dry TC plates



WTS-P100443 Aerobic Bacteria Test Pack (40)

- 40 x 2ml sterile syringes
- 40 x alcohol disinfecting wipes
- 40 x dilution pots
- 40 x compact dry TC plates

WTS-B100817 IDEXX Colilert 18 (20)

WTS-P101401 Coliform & E.Coli Bacteria Test Pack

- 25 x 100ml sterile thio bags
- 20 x IDEXX Colilert 18



WTS-B100023 IDEXX Entrolert DW (20)

WTS-P100445 Enterococci Bacteria Test Pack

- 25 x 100ml sterile thio bags
- 20 x IDEXX Entrolert DW



WTS-B100022 IDEXX Pseudalert (20)

WTS-P100446 Pseudomonas Bacteria Test Pack

- 25 x 100ml sterile thio bags
- 20 x IDEXX Pseudalert



WTS-B100085 Colitag (25)

WTS-P100444 Colitag Bacteria Test Pack
25 x 100ml sterile thio bags
25 x Colitag



WTS-B100840 Readycult® Coliforms 100 (20)

WTS-P101648 Readycult® Coliforms 100 Bacteria Test Pack
25 x 100ml sterile thio bags
20 x Enterococci 100



WTS-B100720 Readycult® Enterococci 100 (20)

WTS-P101649 Readycult® Enterococci 100 Bacteria Test Pack
25 x 100ml sterile thio bags
20 x Enterococci 100

WTS-P101367 IDEXX SimPlate Test Pack
20 x IDEXX SimPlates
2 x 100ml media bottles



WTS-B100099 Legionella Biofilm Kit (5)(Stored Water)

- Test time < 35 Minutes
- Tests Included 5
- Sensitivity 200CFU/Swabbed Area



WTS-B100100 Legionella Industrial Kit (5)(Distribution System)

- Test time < 35 minutes
- Tests included 5
- Sensitivity 100CFU/L



WTS-B100050 Digital Touch Screen Incubator

- External controls
- Clear viewing lid
- Operating temperature: 35 °C - 45 °C
- Dimensions: 215 x 161 x 170 mm
- Power: 100/240v
- Weight: 4 kg



WTS-B100488_UV Lamp (battery operated)



- WTS-B100087 100ml sterile thio bags (25)
- WTS-B100086 100ml sterile thio bags (100)
- WTS-B100364_100ml sterile bottles (25)
- WTS-B100269 100ml sterile bottles (100)

